

Intravet Import Antech Results v.212 and higher Antech Lab Interface

Antech Lab Interface:

- Results can now be downloaded from the Antech website.
- Contact Antech to create your account with username and password.

HOW – TO

- Open website: www.zoasis.com
- Click on Veterinarians link
- Enter Hospital Account #
- Enter Zip Code
- Enter Username
- Enter Password
- Click on “New lab results” link
- Check boxes for all results to be imported
- Click on “Actions” drop-down arrow
- Select “Import”
- Internet Explorer will open a new window of what looks like a lot of programming code- Click on “File” and then “Save As.”
- Rename the file and Save it as an .XML type. (It will default to XML.)

To load results into Intravet:

- In Intravet, select “File” and then “Lab Interface ..Antech”
- In the main window, Click “Open”
- Select the downloaded file
- Click “Ok”
- This will bring all test results into the Antech Lab Interface window. The results will be automatically attached to the right client if there is an exact match in your IntraVet database. Otherwise, you can manually select the correct account.
- Double click on the test you would like to post to pet history—assign a procedure code to the results if you want to send it to invoicing. Otherwise, when you post it, it will go directly to the patient’s history.